

# THD SURGY Mini LightScope



The THD SURGY MINI LED anoscope is a **DISPOSABLE** medical device deemed for investigation of the anal canal.

The THD SURGY MINI LED anoscope combines in a single device a light anoscope and a powerful, fully integrated, LED light source, in order to have a clear and detailed vision of the inspection field. The THD SURGY MINI LED anoscope is ready-to-use and does not need for other tools.

**Bonus:**

- It does not need for an external light source.
- It does not need for a set-up prior the use.
- It avoids the risk of cross-contamination

**SPECIFICATIONS**

REF	Denomination	Width	Length	Window
800108	THD Surgy Mini LED	24 mm	122 mm	85x26 mm



# THD SURGY Mini LightScope

## Specifications

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- Self-light, with an integrated light source
- Ready-to-use
- Light (40g)
- Class I CE device
- Individually packaged
- Latex Free



Ready-to-use integrated light source  
No needs for other devices  
Wide lighting

Fully disposable  
Erase the risk of cross contamination



THD S.p.A. - via Industria, 1 - 42015 Correggio (RE) - Italy  
Tel: 0522634311 Fax: 0522693552 email: [sales@thdlab.com](mailto:sales@thdlab.com)

[www.thdlab.com](http://www.thdlab.com)