



SPON ANAL



Use

SPON ANAL is a compressed surgical haemostatic sponge, manufactured from highly purified neutral gelatine material, for use in ano-rectal surgery.

General Description

SPON ANAL is used during ano-rectal surgical procedures, when the standard haemostatic procedures can not be used or are difficult to apply, or in case that other absorbable materials are inopportune.

SPON ANAL absorb approximately 40-50 times its weight of water and adheres easily to the bleeding site. The uniform porosity of THD SPON ANAL guarantees a favourable haemostasis.



Technical Specifications	
Diameter	30 mm
Length	80 mm

Packaging		
P/N	Qty	Denomination
800050	20	Spon Anal