




THD LIGHTSCOPE RECTOSCOPE LED



Use

The rectoscope THD LightScope LED is a DISPOSABLE medical device, indicated for analysis of ano-rectal canal and the distal region of the sigmoid.

General Description

The rectoscope THD LightScope LED combines into a single disposable device a hand lens and an integrated

powerful LED light source, in order to have a clear and detailed vision of inspection of the area.

The THD Light-Scope rectoscope is ready for use and is also available with the blowing pump included in the sales package.

Features

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- Self-light, with an integrated light source
- 4x magnifier lens
- Two length available 200 mm and 250 mm (7.87 in and 9.84 in)
- Available also with disposable insufflator
- Ready to use
- Light (55 g/0.12 pound)
- CE Class I
- Individually packaged
- Latex Free



Technical Data

P/N. 800097
 P/N. 800098
 P/N. 800099
 P/N. 800100

EN

Advantages

- It doesn't require an additional light source
- It doesn't require a preparation before its use
- It delete the risk of cross-contamination

| Dimensional Features | | |
|--------------------------|------------------|------------------|
| THD LIGHTSCOPE RECTO LED | | |
| | | |
| | Mini | Maxi |
| Diameter (B) | 23 mm (0.90 in) | 23 mm (0.90 in) |
| Length (A) | 200 mm (7.87 in) | 250 mm (9.84 in) |

| THD LightScope Recto LED Sterile | | |
|----------------------------------|-----|--------------------------------------|
| Code | Qty | Denomination |
| 800097 - 10 | 10 | THD LightScope Recto Maxi LED |
| 800098 - 10 | 10 | THD LightScope Recto Mini LED |
| 800099 - 10 | 10 | THD LightScope Recto Maxi LED & bulb |
| 800100 - 10 | 10 | THD LightScope Recto Mini LED & bulb |