




THD LIGHTSCOPE RECTOSCOPE



Use

The rectoscope THD LightScope is a DISPOSABLE medical device, indicated for analysis of ano-rectal canal and the distal region of the sigmoid.

General Description

The rectoscope THD LightScope integrates in the cap an hand lens. The handle has the housing for use with

the light source THD LED2000 and THD REVOLUTION (light source function).

The THD Light-Scope rectoscope is ready-to-use and it is also available with the blowing pump included in the sales package. THD Light-Scope Recto is available with two lengths (200mm and 250mm).

Features

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- 4x magnifier lens
- Two length available 200 mm and 250 mm (7.87 in and 9.84 in)
- Available also with disposable insufflator
- Ready to use
- Light (55 g/0.12 pound)
- CE Class Is
- Individually packaged
- Latex Free



Technical Data

P/N. 800080
P/N. 800081
P/N. 800094
P/N. 800095

EN

Advantages

- It doesn't require a preparation before its use
- It delete the risk of cross-contamination

Dimensional Features		
THD LIGHTSCOPE RECTO		
	Mini	Maxi
Diameter (B)	23 mm (0.90 in)	23 mm (0.90 in)
Length (A)	200 mm (7.87 in)	250 mm (9.84 in)

THD LightScope Recto		
Code	Qty	Denomination
800080 - 10	10	THD LightScope Recto Mini & bulb
800081 - 10	10	THD LightScope Recto Maxi & bulb
800094 - 10	10	THD LightScope Recto Maxi
800095 - 10	10	THD LightScope Recto Mini