



# THD LIGHTSCOPE DOUBLE BULB INSUFFLATOR



### Use

The double bulb insufflator is a device that can be used with THD LightScope rectoscopes line.

It connects quickly to the spout integrated in the THD rectoscope cap.

Usable with devices with spout diameter up to about 6 mm (0.236 in).

The double bulb allows an optimal pressure to be maintained inside the rectoscope, for the best insertion of the device.

The THD LightScope insufflator is reusable

### Features

- Made of soft PVC
- Does not contain phthalates
- Does not contain substances included in the list SVHC
- Does not contain polycyclic aromatic hydrocarbons
- Latex Free
- Confezionato singolarmente
- CE Class I

### Advantages

- Compatible with THD LightScope Rectoscopes
- The lung guarantees pressure maintenance and facilitates the insertion of the device

THD LightScope Double Bulb Insufflator		
Code	Qty	Denomination
800111 - 1	1	THD LightScope Double Bulb insufflator
800111 - 50	50	THD LightScope Double Bulb insufflator